

CATALOG

BEIJING L&Z MEDICAL

Promote scientific innovation

Stand up to future challenges

Become a global leading medical company



||| Company Profile

- Founded in 2001, with over 300 employees, including overseas employees
- Headquarters located in Beijing
- One of the leaders in Chinese medical field
- Provides quality products and best solutions to patients and society
- Full range of innovative Enteral and Parenteral feeding devices
- 33% market share in Parenteral and Enteral feeding devices market in China
- Own a number of patents for invention and for utility model
- Industrial layout: China & US factories, China & US R&D center

Our vision

Actively apply scientific innovation, calmly meet future challenges, strive to become a leading global medical device enterprise

Our mission

Provide innovative medical solutions for patients and society

Core values

Care for life, scientific innovation, keep improving



Beijing L&Z Medical Technology Development Co., Ltd


Telephone: 010-87952982 010-87952983

Fax: 010-87952199

E-mail: info@lingze.com

Address: No. 8, M2-5 Block, Xinggu Industrial Developing Zone,
Pinggu District 101200 Beijing China


||| Disposable Enteral Feeding Sets



Bag Pump set (use with pump)

Model: 300ml, 500ml, 600ml, 1000ml, 1200ml, 1500ml


||| Disposable Enteral Feeding Sets




Bag Gravity set

Model: 300ml, 500ml, 600ml, 1000ml, 1200ml, 1500ml

||| Disposable Enteral Feeding Sets



||| Disposable Enteral Feeding Sets



Spike Gravity set

Model: Vented spike, Non-vented spike
Non-vented ENPlus spike, Universal ENPlus spike

||| Disposable Enteral Feeding Sets



Double bag set

Model: 1000/1200ml feeding bag+1000/1200ml flushing bag
1000/1200ml flushing bag+spike/screw spike

||| Nasogastric tubes

Y port option ←



Nasogastric tube

PVC - CH12/14/16, 120cm
PUR - CH8/10/12/14/16, 110/130/150cm

Enteral Feeding Pump



Characteristic			
Infusion rate range	0-400ml/h(increment 1ml/h), with error less than $\pm 10\%$	Temperature	+10°C ~ +30°C
Task setting range	0-3000ml (increment 1ml)	Relative humidity	30% ~ 75%
Dosing range	0-3000ml (increment 1ml)	Air pressure	700hpa ~ 1060hpa
Interval setting range	1h-24h (increment 1h)	Language	English/Chinese
Total infusion volume	0-9999ml (minimum 1ml displayable)	Power supply	100V-240V, 50/60HZ
Clock	Built-in 24h time clocking	Memory function	24h
Alarm volume	Low, middle and high	Battery cruising capacity	Non-stop 24h under the pump operating at an infusion rate of 25ml/h
Screen	LCD	Economic mode	Screen off mode
Battery	4 chargeable lithium cells, DC8.4V	Service life	5 years

Enteral feeding pump

With or W/O heating system

||| Disposable infusion bag for parenteral nutrition



TPN bag

Different configurations could be chosen



||| Peripherally inserted central catheter

4F Single Lumen

5F Dual Lumen

PICC

Made in US



Enteral and Parenteral
nutrition feeding
specialist



Beijing L&Z Medical Technology
Development Co., Ltd
Telephone: 010-87952982 010-87952983
Fax: 010-87952199
E-mail: info@lingze.com
Address: No. 8, M2-5 Block, Xinggu
Industrial Developing Zone, Pinggu
District 101200 Beijing China

BEIJING L&Z MEDICAL TECHNOLOGY DEVELOPMENT CO., LTD