




# GLUTTON LIFE BIOPSY FORCEPS



SAMPLE RETRIEVER  
INCLUDED

SECURITY & ERGONOMY



**Description:**

- PTFE coated biopsy forceps,
- Graduated catheter to easily identify the biopsy forceps progression whitin the working channel.
- All forceps are delivered with a sample retriever to reduce blood exposure risk.

Primary packaging: individual pouch  
Secondary packaging: cardboard box

References	Biopsy forceps	Forceps		Catheter		Minimum delivery quantity
		Jaws	Spike	Length (cm)	Diameter (mm)	
CJ-PPT-18-110-A	Tracheo-bronchial	= standard	Yes	110	1.8	50
CJ-PPT-18-160	Pediatric	= standard	No	160	1.8	20
CJ-PPT-18-160-A	Pediatric	= standard	Yes	160	1.8	20
CJ-PAT-23-160	Gastric	= standard	No	160	2.3	50
CJ-PAT-23-160-A	Gastric	= standard	Yes	160	2.3	50
CJ-PATC-23-160	Gastric	≡ alligator	No	160	2.3	50
CJ-PATC-23-160-A	Gastric	≡ alligator	Yes	160	2.3	50
CJ-PAT-23-230	Gastric and colonic	= standard	No	230	2.3	50
CJ-PAT-23-230-A	Gastric and colonic	= standard	Yes	230	2.3	50
CJ-PATC-23-230	Gastric and colonic	≡ alligator	No	230	2.3	50
CJ-PATC-23-230-A	Gastric and colonic	≡ alligator	Yes	230	2.3	50