


DILATION & ACHALASIA BALLOONS

SINGLE DIAMETER



BACK TO SUMMARY



Primary packaging: Pouch

Secondary packaging: cardboard box

Description:

- Available for: oesophagus, cardia, pylorus, colon, biliary duct,
- One size diameter,
- Dilation balloons are TTS compatible,
- Achalasia balloons are OTW procedure compatible.
- Dilation Balloons are compatible with inflation devices NLS/ID30A, VTL/ID30B and VTL/ID-60 while achalasia balloons are not.

References	Balloon		Catheter		Pressure (ATM)		Comments	Minimum delivery quantity
	Diameter (mm)	Length (mm)	Diameter (mm/fr)	Length (cm)	Recommended	Burst		
Dilation balloons								
CJ-JHY-BD-06-40-250	6	40	2.4 / 7.5	250	8	10	Working channel diameter: ≥ 2,8 mm (except CJ-JHY-BD-20-80-250 : ≥ 3.2 mm) Guide wire: ≤ 0,035"	3
CJ-JHY-BD-08-40-250	8	40	2.4 / 7.5	250	8	10		
CJ-JHY-BD-10-40-250	10	40	2.4 / 7.5	250	8	10		
CJ-JHY-BD-12-60-250	12	60	2.4 / 7.5	250	5	6		
CJ-JHY-BD-14-60-250	14	60	2.4 / 7.5	250	5	6		
CJ-JHY-BD-16-60-250	16	60	2.4 / 7.5	250	5	6		
CJ-JHY-BD-16-80-250	16	80	2.4 / 7.5	250	5	6		
CJ-JHY-BD-18-80-250	18	80	2.4 / 7.5	250	4	5		
CJ-JHY-BD-20-80-250	20	80	2.4 / 7.5	250	4	5		
Achalasia balloons								
CJ-JHY-BD-30-80-90	30	80	3.4 / 10	90	2	3	OTW procedure Guide wire: ≤ 0,035"	1
CJ-JHY-BD-35-80-90	35	80	3.4 / 10	90	2	3		
CJ-JHY-BD-40-80-90	40	80	3.4 / 10	90	2	3		