



Disposable Introducer Kit

Quarter-Turn-Screw Dilator Locking
Provide safe and fast system insertion
and easy dilator withdrawal



Smooth Transition
Minimized trauma to the
patients ensures safety

Product Information

Features

- Well balanced pushability and flexibility
- Gap/kink resistance
- Smooth transition from dilator to sheath
- Special designed dilator locking

Benefits

- Easy entry and smoother navigation
- Good hemostatic performance
- Minimized trauma to patients

Ranges

- Femoral Sheath: 4-9F
- Radial Sheath: 4-7F

Kit Components

- Introducer Sheath (1)
- Dilator (1)
- Stainless Steel Guide Wire (1)
- Introducer Needle (1)



Ordering Information

Catalogue No.	Sheath Size	Sheath Length	Radial Sheath	Femoral Sheath
MMQGB04110	4F	11cm	✓	✓
MMQGB05110	5F	11cm	✓	✓
MMQGB06110	6F	11cm	✓	✓
MMQGB07110	7F	11cm	✓	✓
MMQGB08110	8F	11cm		✓
MMQGB09110	9F	11cm		✓
MMQGB04160	4F	16cm	✓	✓
MMQGB05160	5F	16cm	✓	✓
MMQGB06160	6F	16cm	✓	✓
MMQGB07160	7F	16cm	✓	✓
MMQGB08160	8F	16cm		✓
MMQGB09160	9F	16cm		✓
MMQGB04240	4F	24cm	✓	✓
MMQGB05240	5F	24cm	✓	✓
MMQGB06240	6F	24cm	✓	✓
MMQGB07240	7F	24cm	✓	✓
MMQGB08240	8F	24cm		✓
MMQGB09240	9F	24cm		✓

Please contact the regional sales for more information about this product in your country.