


NOVAJECT INJECTION NEEDLE



BACK TO SUMMARY



Primary packaging: Pouch
Secondary packaging: cardboard box

Description:

- Triple bevel needle for an easy penetration,
- Ergonomic handle to control needle movements,
- Transparent catheter for liquid flow visualization.

References	Description	Catheter		Needle		Handle color	Minimum delivery quantity
		Diameter (mm/fr)	Length (cm)	Diameter (G / mm)	Length (mm)		
CJ-A3-23-230-21*	Novaject	2.3	230	21 / 0.8	3	●	10
CJ-A3-23-230-23*		2.3	230	23 / 0.6	3	●	
CJ-A3-23-230-25*		2.3	230	25 / 0.5	3	●	
CJ-A-18-165-25		1.8	165	25 / 0.5	4	●	
CJ-A-23-160-19		2.3	160	19 / 1.1	4	●	
CJ-A-23-230-19		2.3	230	19 / 1.1	4	●	
CJ-A-23-230-21		2.3	230	21 / 0.8	4	●	
CJ-A-23-230-23		2.3	230	23 / 0.6	4	●	
CJ-A-23-230-25		2.3	230	25 / 0.5	4	●	
CJ-A5-23-230-21*		2.3	230	21 / 0.8	5	●	
CJ-A5-23-230-23		2.3	230	23 / 0.6	5	●	
CJ-A5-23-230-25		2.3	230	25 / 0.5	5	●	

* Special request, delivery time on demand