

Angiopower™

Inflation Device



Product Information

Inflation Device:

Features

- Analog/fluorescent display
- Symmetrical handling
- Ergonomic design
- Barrier screw pivot
- Smooth fine adjustment
- Rapid trigger release

Benefits

- Visual pressure control
- For right and left hand users
- Easy piston grip
- Rapid inflation
- Effective pressure transmission
- Rapid/secure deflation

Ranges

- Maximum pressure: 30atm
- Barrel volume: 20mL
- Stopcock: 3-way 500psi
- Tubing: length 33cm 1200psi



Technical Parameters

Display	Analog/Fluorescent
Pressure	30atm (440psi)
Barrel Volume	20mL
Stopcock	3-way 500psi
Extension Tubing	33cm, 1200psi, male rotating Luer connector

Ordering Information

Catalogue No.	Inflation Range	Stopcock
LP-P-30	30 atm	/
LP-P-30S	30 atm	Included



Please contact the regional sales for more information about this product in your country.

37 Chaoqian Road, Changping District, Beijing, China 102200, +86 10 80120666, international@lepu-medical.com, en.lepumedical.com

Inflation Devices Kit



Product Information

Inflation Device:

Features

- Analog/fluorescent display
- Symmetrical handling
- Ergonomic design
- Barrier screw pivot
- Smooth fine adjustment
- Rapid trigger release

Benefits

- Visual pressure control
- For right and left hand users
- Easy piston grip
- Rapid inflation
- Effective pressure transmission
- Rapid/secure deflation

Ranges

- Maximum pressure: 30atm
- Barrel volume: 20mL
- Stopcock: 3-way 500psi
- Tubing: length 33cm 1200psi

Y-Hemostasis Valve Set:

Ranges

- Rotating / Push-Pull 250psi
- Guide Wire Introducer
- Guide Wire Torquer

Ordering Information

Catalogue No.	Inflation Device Model		Y-Hemostasis Valve Set Model			
	Type	Inflation Range	Type	Mode	Contains	Extension Tubing
LP-P-30-YNP	LP-P-30	30atm	YNP	Push-Pull	Whole set	/
LP-P-30-YNP20	LP-P-30	30atm	YNP20	Push-Pull	Whole set	20cm
LP-P-30-YNP25	LP-P-30	30atm	YNP25	Push-Pull	Whole set	25cm
LP-P-30-YNP40	LP-P-30	30atm	YNP40	Push-Pull	Whole set	40cm
LP-P-30-YNP50	LP-P-30	30atm	YNP50	Push-Pull	Whole set	50cm
LP-P-30-YNPH	LP-P-30	30atm	YNPH	Push-Pull	Y-valve only	/
LP-P-30-YNP20H	LP-P-30	30atm	YNP20H	Push-Pull	Y-valve only	20cm
LP-P-30-YNP25H	LP-P-30	30atm	YNP25H	Push-Pull	Y-valve only	25cm
LP-P-30-YNP40H	LP-P-30	30atm	YNP40H	Push-Pull	Y-valve only	40cm
LP-P-30-YNP50H	LP-P-30	30atm	YNP50H	Push-Pull	Y-valve only	50cm
LP-P-30-YNR	LP-P-30	30atm	YNR	Rotating	Whole set	/
LP-P-30-YNRH	LP-P-30	30atm	YNRH	Rotating	Y-valve only	/
LP-P-30S-YNP	LP-P-30S	30atm	YNP	Push-Pull	Whole set	/
LP-P-30S-YNP20	LP-P-30S	30atm	YNP20	Push-Pull	Whole set	20cm
LP-P-30S-YNP25	LP-P-30S	30atm	YNP25	Push-Pull	Whole set	25cm
LP-P-30S-YNP40	LP-P-30S	30atm	YNP40	Push-Pull	Whole set	40cm
LP-P-30S-YNP50	LP-P-30S	30atm	YNP50	Push-Pull	Whole set	50cm
LP-P-30S-YNPH	LP-P-30S	30atm	YNPH	Push-Pull	Y-valve only	/
LP-P-30S-YNP20H	LP-P-30S	30atm	YNP20H	Push-Pull	Y-valve only	20cm
LP-P-30S-YNP25H	LP-P-30S	30atm	YNP25H	Push-Pull	Y-valve only	25cm
LP-P-30S-YNP40H	LP-P-30S	30atm	YNP40H	Push-Pull	Y-valve only	40cm
LP-P-30S-YNP50H	LP-P-30S	30atm	YNP50H	Push-Pull	Y-valve only	50cm
LP-P-30S-YNR	LP-P-30S	30atm	YNR	Rotating	Whole set	/
LP-P-30S-YNRH	LP-P-30S	30atm	YNRH	Rotating	Y-valve only	/