



Shoocin™

Introducer Kit



Quarter-Turn-Screw Dilator Locking

Provide safe and fast system insertion and easy dilator withdrawal



Smooth Transition

Minimized trauma to the patients ensures safety

Product Information

Features

- Well balanced pushability and flexibility
- Gap/kink resistance
- Smooth transition from dilator to sheath
- Special designed dilator locking

Benefits

- Easy entry and smoother navigation
- Good hemostatic performance
- Minimized trauma to patients

Ranges

- Femoral Sheath: 4-9F
- Radial Sheath: 4-7F
- Long Sheath: 6-8F

Full Kit Components

- Introducer Sheath (1)
(radial/femoral/long sheath)
- Stainless Steel Guide Wire (1)
- Introducer Needle (1)
- Scalpel (1)
- Syringe (1)

Ordering Information*

Femoral Sheath

Catalogue No. (Stainless Steel Guide Wire)	Sheath Size	Sheath Length	Guide Wire Diameter	Guide Wire Length	Introducer Needle	Scalpel	Syringe
FS040525-SS45-IN20-K-S	4F	5cm	0.025"	45cm	20G	1	1
FS041125-SS45-IN20-K-S	4F	11cm	0.025"	45cm	20G	1	1
FS041135-SJ45-IN18-K-S	4F	11cm	0.035"	45cm	18G	1	1
FS050535-SJ45-IN18-K-S	5F	5cm	0.035"	45cm	18G	1	1
FS050538-SJ45-IN18-K-S	5F	5cm	0.038"	45cm	18G	1	1
FS051135-SJ45-IN18-K-S	5F	11cm	0.035"	45cm	18G	1	1
FS051138-SJ45-IN18-K-S	5F	11cm	0.038"	45cm	18G	1	1
FS060535-SJ45-IN18-K-S	6F	5cm	0.035"	45cm	18G	1	1
FS060538-SJ45-IN18-K-S	6F	5cm	0.038"	45cm	18G	1	1
FS061135-SJ45-IN18-K-S	6F	11cm	0.035"	45cm	18G	1	1
FS061138-SJ45-IN18-K-S	6F	11cm	0.038"	45cm	18G	1	1
FS070535-SJ45-IN18-K-S	7F	5cm	0.035"	45cm	18G	1	1
FS070538-SJ45-IN18-K-S	7F	5cm	0.038"	45cm	18G	1	1
FS071135-SJ45-IN18-K-S	7F	11cm	0.035"	45cm	18G	1	1
FS071138-SJ45-IN18-K-S	7F	11cm	0.038"	45cm	18G	1	1
FS080535-SJ45-IN18-K-S	8F	5cm	0.035"	45cm	18G	1	1
FS080538-SJ45-IN18-K-S	8F	5cm	0.038"	45cm	18G	1	1
FS081135-SJ45-IN18-K-S	8F	11cm	0.035"	45cm	18G	1	1
FS081138-SJ45-IN18-K-S	8F	11cm	0.038"	45cm	18G	1	1
FS090538-SJ45-IN18-K-S	9F	5cm	0.038"	45cm	18G	1	1
FS091138-SJ45-IN18-K-S	9F	11cm	0.038"	45cm	18G	1	1

Ordering Information*

Radial Sheath

Catalogue No. (Stainless Steel Guide Wire)	Sheath Size	Sheath Length	Guide Wire Diameter	Guide Wire Length	Introducer Needle	Scalpel	Syringe
RS040518-SS45-IN20-K-S	4F	5cm	0.018"	45cm	20G	1	1
RS040521-SS45-IN20-K-S	4F	5cm	0.021"	45cm	20G	1	1
RS040525-SS45-IN20-K-S	4F	5cm	0.025"	45cm	20G	1	1
RS040718-SS45-IN20-K-S	4F	7cm	0.018"	45cm	20G	1	1
RS040721-SS45-IN20-K-S	4F	7cm	0.021"	45cm	20G	1	1
RS040725-SS45-IN20-K-S	4F	7cm	0.025"	45cm	20G	1	1
RS050518-SS45-IN20-K-S	5F	5cm	0.018"	45cm	20G	1	1
RS050521-SS45-IN20-K-S	5F	5cm	0.021"	45cm	20G	1	1
RS050525-SS45-IN20-K-S	5F	5cm	0.025"	45cm	20G	1	1
RS050718-SS45-IN20-K-S	5F	7cm	0.018"	45cm	20G	1	1
RS050721-SS45-IN20-K-S	5F	7cm	0.021"	45cm	20G	1	1
RS050725-SS45-IN20-K-S	5F	7cm	0.025"	45cm	20G	1	1
RS060518-SS45-IN20-K-S	6F	5cm	0.018"	45cm	20G	1	1
RS060521-SS45-IN20-K-S	6F	5cm	0.021"	45cm	20G	1	1
RS060525-SS45-IN20-K-S	6F	5cm	0.025"	45cm	20G	1	1
RS060718-SS45-IN20-K-S	6F	7cm	0.018"	45cm	20G	1	1
RS060721-SS45-IN20-K-S	6F	7cm	0.021"	45cm	20G	1	1
RS060725-SS45-IN20-K-S	6F	7cm	0.025"	45cm	20G	1	1
RS041121-SS45-IN20-K-S	4F	11cm	0.021"	45cm	20G	1	1
RS042421-SS70-IN20-K-S	4F	24cm	0.021"	70cm	20G	1	1
RS051121-SS45-IN20-K-S	5F	11cm	0.021"	45cm	20G	1	1
RS052421-SS70-IN20-K-S	5F	24cm	0.021"	70cm	20G	1	1
RS061121-SS45-IN20-K-S	6F	11cm	0.021"	45cm	20G	1	1
RS062421-SS70-IN20-K-S	6F	24cm	0.021"	70cm	20G	1	1
RS070521-SS45-IN20-K-S	7F	5cm	0.021"	45cm	20G	1	1
RS070521-SS70-IN20-K-S	7F	5cm	0.021"	70cm	20G	1	1
RS071121-SS45-IN20-K-S	7F	11cm	0.021"	45cm	20G	1	1
RS072421-SS70-IN20-K-S	7F	24cm	0.021"	70cm	20G	1	1
RS070718-SS45-IN20-K-S	7F	7cm	0.018"	45cm	20G	1	1
RS070718-SS70-IN20-K-S	7F	7cm	0.018"	70cm	20G	1	1
RS071118-SS45-IN20-K-S	7F	11cm	0.018"	45cm	20G	1	1
RS071118-SS70-IN20-K-S	7F	11cm	0.018"	70cm	20G	1	1

Long Sheath

Catalogue No. (Stainless Steel Guide Wire)	Sheath Size	Sheath Length	Guide Wire Diameter	Guide Wire Length	Introducer Needle	Scalpel	Syringe
LS061638-SJ80-IN18-K-S	6F	16cm	0.038"	80cm	18G	1	1
LS062338-SJ80-IN18-K-S	6F	23cm	0.038"	80cm	18G	1	1
LS071638-SJ80-IN18-K-S	7F	16cm	0.038"	80cm	18G	1	1
LS072338-SJ80-IN18-K-S	7F	23cm	0.038"	80cm	18G	1	1
LS081638-SJ80-IN18-K-S	8F	16cm	0.038"	80cm	18G	1	1
LS082338-SJ80-IN18-K-S	8F	23cm	0.038"	80cm	18G	1	1

* The Yellow Marks indicate the components which can be customized as following:

- Remove the every single component includes guide wire, introducer needle (IN), scalpel (K) and syringe (S).
- Replace the stainless steel guide wire (straight SS/J tip SJ) into hydrophilic guide wire (straight HS/135° angled HW) or nitinol guide wire (straight NS/J tip NJ).
- Replace the introducer needle (IN) into intravascular catheter (IC). Intravascular catheter is mandatory to be a set with hydrophilic guide wire.

(e.g. Radial sheath/6F/7cm/0.021" straight stainless steel guide wire/45cm/20G introducer needle/without scalpel and syringe: RS060721-SS45-IN20)

Please contact the sales for more information about the selections of customized kit.