

# Angioway™

Manifolds **CE**<sub>0197</sub>




## Ordering Information

Catalogue No.	Pressure	Side Port	Valve Type	Extension Tubing
LPJD502N	500psi	2	On	/
LPJD502F	500psi	2	Off	/
LPJD503N	500psi	3	On	/
LPJD503F	500psi	3	Off	/
LPJD503NX	500psi	3	On	Included
LPJD503FX	500psi	3	Off	Included

# Angiospring™

Angiography Control Syringe **CE**<sub>0197</sub>



## Ordering Information

Catalogue No.	Connector Type	Nominal Capacity
LPHB112	Rotating	12 ml

## Product Information


### Features

- Side port: 2/3
- Pressure: 500psi
- Valve type: On / Off
- Extension tubing: optional

### Features

- Solid and smooth plunger
- Clear scale
- Latex-free

# Contrast Medium Injection Manifold Kits



## Product Information

### Manifold Kit Includes

- Manifold (1)
- Angiography Control Syringe (1)
- High Pressure Extension Tubing (1)
- Infusion Set (2)

## Ordering Information

Catalogue No.	Manifold Model	Control Syringe Model	High Pressure Extension Tubing Model	Infusion Set Model	
MK-005	LPJD253N	LPHB112	LPJG30	IS-G-V	IS-G-V3
MK-006	LPJD253N	LPHB112	LPJG60	IS-G-V	IS-G-V3
MK-007	LPJD253N	LPHB112	LPJG90	IS-G-V	IS-G-V3
MK-008	LPJD253N	LPHB112	LPJG120	IS-G-V	IS-G-V3
MK-017	LPJD253F	LPHB112	LPJG30	IS-G-V	IS-G-V3
MK-018	LPJD253F	LPHB112	LPJG60	IS-G-V	IS-G-V3
MK-019	LPJD253F	LPHB112	LPJG90	IS-G-V	IS-G-V3
MK-020	LPJD253F	LPHB112	LPJG120	IS-G-V	IS-G-V3