




NOVA CLIP R² HEMOCLIP



Primary packaging: individual pouch
Secondary packaging: cardboard box

Description:

- Innovative handle enabling a 360° rotation with one hand,
- Repositionnable,
- Optimal flexibility of the distal part,
- Ready for use,
- MRI conditional.

References	Clip			Catheter		Handle color	Minimum delivery quantity
	Maximal opening (mm)	Jaws length (mm)	Jaws angle	Diameter (mm)	Length (cm)		
CJ-HC09-25-230	9	5	135°	2.5	230	●	10
CJ-HC11-25-230	11	6	135°	2.5	230	●	10
CJ-HC11-25-160*	11	6	135°	2.5	160	●	10
CJ-HC11-25-260	11	6	135°	2.5	260	●	10
CJ-HC13-25-230	13	7	135°	2.5	230	●	10
CJ-HC1690-25-230	16	9	90°	2.5	230	●	10

Minimum working channel diameter : 2.8 mm

* Special request, delivery time on demand